

Upper Extremity Superficial Venous Mapping Image Summary

******These images should be taken with a split screen showing the vein compression on the right and the AP diameter measurement on the left side of screen.

<u>Image</u>	<u>Measurement</u>
1. Patient information screen	
2. Trans dual screen Cephalic prox upper arm	AP diameter
3. Trans dual screen Cephalic mid upper arm	AP diameter
4. Trans dual screen Cephalic dist upper arm	AP diameter
5. Trans dual screen Cephalic prox forearm	AP diameter
6. Trans dual screen Cephalic mid forearm	AP diameter
7. Trans dual screen Cephalic dist forearm	AP diameter
8. Trans dual screen Basilic prox upper arm	AP diameter
9. Trans dual screen Basilic mid upper arm	AP diameter
10. Trans dual screen Basilic dist upper arm	AP diameter
11. Trans dual screen Basilic prox forearm	AP diameter
12. Trans dual screen Basilic mid forearm	AP diameter
13. Trans dual screen Basilic dist forearm	AP diameter
14. If veins appear small (<2.5mm), apply a tourniquet and repeat measurements.	
15. Repeat images in contralateral arm if a bilateral exam is indicated.	

Always document if thrombus or wall thickening is seen.