

Creating Customer User Specific Window Level Settings Diagnostic Desktop

Do NOT modify an existing window level preset. Follow the instructions below to create a NEW window level preset.

1. Click the dropdown arrow on the window level tool to access the Edit button.
2. Click the Edit button.
3. Click the New window level preset button.
4. Type the name for your new window level (W/L) preset with your initials at the end.
5. Set the Window width.
6. Set the Window center.
7. Set the filter (if one is needed).
8. Use the up or down arrow to move your NEW window level (W/L) setting to the preset number that you want it to correspond to.
9. Click OK.

