

Pediatric Renal Ultrasound Image Summary

1. pancreas – trans (head / uncinata)
2. pancreas - trans (body/tail)
3. pancreas – long
4. aorta – long prox 1. AP
5. aorta – long mid 2. AP
6. aorta – long distal 3. AP
7. IVC long
8. IVC long (color)
9. right kidney/liver – long
10. right kidney – long 8. Length
11. right kidney - long
12. right kidney - with color
13. right kidney - long medial
14. right kidney - long lateral
15. right kidney - trans superior
16. right kidney - trans mid
17. right kidney - trans mid
18. right kidney – trans inferior

19. left kidney/ spleen – long
20. left kidney – long 18. Length
21. left kidney – long
22. left kidney – with color
23. left kidney – long medial
24. left kidney – long lateral
25. left kidney - trans superior
26. left kidney - trans mid
27. left kidney – trans mid
28. left kidney - trans inferior

27. bladder – long 27 . AP & length
28. bladder – trans 28. AP & transverse

29. Assess the bladder wall with a high frequency linear probe (12MHZ)

30. ureteral jet documentation (color)
31. bladder (post void) – long 31. AP & length
32. bladder (post void) – trans 32. AP & transverse